



CUSTOM PACKAGING + SPECIALTY BAGS

Tailored Packaging Solutions for Your Organization

Cenmed Enterprises, with its 30 years of experience in healthcare, research, and science, understands the critical role packaging plays in your organization's success. They offer a unique advantage: custom-designed bags and boxes that perfectly align with your specific needs and requirements.



By partnering with Cenmed, you'll gain access to custom bags and boxes that not only meet your specific needs but also enhance your organization's efficiency, professionalism, and brand identity.

- **Research Labs:** Custom-designed specimen collection kits with compartments for different types of samples, temperature control features, and clear labeling for easy identification.
- **Hospitals and Clinics:** Secure, tamper-evident medication pouches personalized with patient information and dosage instructions.
- **Medical Equipment Companies:** Sturdy, reusable transport containers for delicate equipment, with cushioning and compartments for secure storage.

Supporting your mission in
Healthcare | Research | Science™

Amber Zip Bags



- **Blocks harmful UV light:** Shield light-sensitive medications during transport.
- **Pharmacy to patient, leak-proof:** Safely route meds with these 2-mil bags.
- **Bulk convenience:** Packed 100 per case for efficient restocking.

Ordering Information

Closure	Color	Image	Material	Width x Length	Thickness	Text
Zip Seal Closure	Amber		LDPE	2W x 3L	2 Mil	
Zip Seal Closure	Amber		LDPE	2.5W x 9L	2 Mil	
Zip Seal Closure	Amber		LDPE	3W x 5L	2 Mil	
Zip Seal Closure	Amber		LDPE	4W x 6L	2 Mil	
Zip Seal Closure	Amber		LDPE	6W x 9L	2 Mil	
Zip Seal Closure	Amber		LDPE	8W x 14L	2 Mil	
Zip Seal Closure	Amber		LDPE	9W x 12L	2 Mil	
Zip Seal Closure	Amber		LDPE	12W x 12L	2 Mil	

CUSTOM PACKAGING + SPECIALTY BAGS

Autoclave Bags



- **Sterilize infectious waste safely:** High-temperature autoclave bags eliminate harmful pathogens.
- **Handles intense heat:** Perfect for biohazardous materials up to 275°F (135°C).
- **Clear sterilization confirmation:** Biohazard symbol and "Autoclaved" indicator ensure successful treatment.

Ordering Information

Closure	Color	Image	Material	Width x Length	Thickness	Text
Twist-Tie Closure	Red	Biohazard	Polyethylene	8W x 12L	2 Mil	Biohazard
Twist-Tie Closure	Red	Biohazard	Polyethylene	12W x 24L	2 Mil	Biohazard
Twist-Tie Closure	Red	Biohazard	Polyethylene	14W x 19L	2 Mil	Biohazard
Twist-Tie Closure	Red	Biohazard	Polyethylene	19W x 23L	2 Mil	Biohazard
Twist-Tie Closure	Red	Biohazard	Polyethylene	24W x 30L	2 Mil	Biohazard
Twist-Tie Closure	Red	Biohazard	Polyethylene	25W x 35L	2 Mil	Biohazard

Chemotherapy Bags




- **Secure transport:** 4-mil thick, reclosable bags safely contain cytotoxic drugs, syringes, and more.
- **Clear communication:** Bilingual biohazard warnings printed in black and yellow for immediate recognition.
- **Convenient storage:** Packed in 100-count cases for easy restocking and organization.

Ordering Information

Closure	Color	Image	Material	Width x Length	Thickness	Text
Zip Seal Closure	Yellow	Chemotherapy	Polyethylene	6W x 9L	2 Mil	Caution
Zip Seal Closure	Yellow	Chemotherapy	Polyethylene	9W x 12L	2 Mil	Caution
Zip Seal Closure	Yellow	Chemotherapy	Polyethylene	12W x 15L	2 Mil	Caution
Zip Seal Closure	Yellow	Chemotherapy	Polyethylene	13W x 19L	2 Mil	Caution

CUSTOM PACKAGING + SPECIALTY BAGS

Specimen Transport Bags



- **Lab efficiency boost:** Streamline workflows and minimize errors with the right specimen transport bag.
- **Double protection:** Leakproof pouch secures specimens, while an open-end pocket keeps paperwork dry.
- **Customize for your lab:** Choose colors, imprints, and sizes to match your needs.

Ordering Information

Closure	Color	Image	Material	Width x Length	Thickness	Text
Zip Seal Closure	Clear	Biohazard	Polyethylene	4W x 6L	1.7 Mil	Biohazard
Zip Seal Closure	Clear	Biohazard	Polyethylene	6W x 6L	1.75 Mil	Biohazard
Zip Seal Closure	Clear	Biohazard	Polyethylene	6W x 9L	1.75 Mil	Biohazard
Zip Seal Closure	Clear	Biohazard	Polyethylene	6W x 9L	2 Mil	Biohazard
Zip Seal Closure	Clear	Biohazard	Polyethylene	7W x 9L	1.75 Mil	Biohazard
Zip Seal Closure	Clear	Biohazard	Polyethylene	7W x 9.5L	2 Mil	Biohazard
Zip Seal Closure	Clear	Biohazard	Polyethylene	8W x 10L	1.75 Mil	Biohazard

Patient Belonging Bags



- **Streamline supply management:** Transport & organize healthcare essentials.
- **Secure your supplies:** Choose from four closure styles - Draw Tape, Drawstring, Rigid Handle, or Tamper Evident.
- **Find your perfect match:** Three colors with custom options available.

Ordering Information

Closure	Color	Image	Material	Width x Length	Thickness	Text
Drawstring Closure	Blue			20W x 20L		Pers. Belongings
Drawstring Closure	White			20W x 20L		Pers. Belongings
Drawstring Closure	Clear			20W x 20L		Pers. Belongings
Drawstring Closure	Clear			20W x 4D x 24L		Pers. Belongings

CUSTOM PACKAGING + SPECIALTY BAGS

Infectious Waste Liners



- **Lab efficiency boost:** Streamline workflows and minimize errors with the right specimen transport bag.
- **Double protection:** Leakproof pouch secures specimens, while an open-end pocket keeps paperwork dry.
- **Customize for your lab:** Choose colors, imprints, and sizes to match your needs.

Ordering Information

Closure	Color	Image		Width x Length	Thickness	Text
X-Seal Closure	Red	Biohazard	Polyethylene	11W x 14L	1.5 Mil	Infectious Waste
X-Seal Closure	Red	Biohazard	Polyethylene	12W x 8D x 30L	1.1 Mil	Infectious Waste
X-Seal Closure	Red	Biohazard	Polyethylene	24W x 23L	1.3 Mil	Infectious Waste
X-Seal Closure	Red	Biohazard	Polyethylene	24W x 24L	11 Micron	Infectious Waste
X-Seal Closure	Red	Biohazard	Polyethylene	24W x 24L	1.5 Mil	Infectious Waste
X-Seal Closure	Red	Biohazard	Polyethylene	24W x 30L	13 Micron	Infectious Waste
X-Seal Closure	Red	Biohazard	Polyethylene	30W x 36L	1.5 Mil	Infectious Waste
X-Seal Closure	Red	Biohazard	Polyethylene	31W x 43L	17 Micron	Infectious Waste
X-Seal Closure	Blue	Biohazard	Polyethylene	31W x 43L	17 Micron	Infectious Waste
X-Seal Closure	Yellow	Biohazard	Polyethylene	31W x 43L	1.2 Mil	Infectious Waste
X-Seal Closure	Red	Biohazard	Polyethylene	31W x 43L	1.5 Mil	Infectious Waste
X-Seal Closure	Clear	Biohazard	Polyethylene	33W x 39L	.06 Mil	Infectious Waste
X-Seal Closure	Red	Biohazard	Polyethylene	36W x 48L	17 Micron	Infectious Waste
X-Seal Closure	Red	Biohazard	Polyethylene	37W x 50L	1.5 Mil	Infectious Waste
X-Seal Closure	Red	Biohazard	Polyethylene	40W x 46L	1.3 Mil	Infectious Waste
X-Seal Closure	Red	Biohazard	Polyethylene	40W x 46L	1.5 Mil	Infectious Waste

CUSTOM PACKAGING + SPECIALTY BAGS

Drug Tray Security Bags



- **Seal-safe guaranteed:** Rest assured knowing this product uses polyethylene approved by both the FDA and USDA.
- **Tamper-proof security:** The secure, tamper-evident tape closure ensures all components arrive untouched and protected.
- **Easy-open design:** A convenient center perforation on the clear 2-mil bag allows for quick and hassle-free access in emergencies.
- **Customizable color options:** Choose between blue or white peel-back paper tabs for a touch of personalized branding or compliance.

Ordering Information

Closure	Color	Image	Material	Width x Length	Thickness	Text
Tamper Evident Tape	Clear		Polyethylene	14W x 22L	2 Mil	
Tamper Evident Tape	Clear		Polyethylene	20W x 29L	2 Mil	

Gas Sterilization Bags



- **Grab & Go Convenience:** Each bag is securely mounted on a stapled header for easy dispensing and on-the-go access. No more fumbling or searching for individual bags.
- **Instantly Visible:** The distinctive green tint ensures quick identification, saving you precious time in fast-paced environments.
- **Bulk Savings & Efficiency:** Each pack comes with 100 pre-mounted bags, minimizing restocking needs and maximizing efficiency.

Ordering Information

Closure	Color	Image	Material	Width x Length	Thickness	Text
	Green			12W x 15L	2 Mil	

CUSTOM PACKAGING + SPECIALTY BAGS

Green IV Covers



- **UV-Blocking Green Tint:** Safeguards light-sensitive contents from harmful ultraviolet rays.
- **Convenient IV Access:** Slit design at the sealed end allows for easy insertion of IV bags.
- **Bulk Packaging:** Inner packs of 100 bags ensure you have what you need, when you need it.

Ordering Information

Closure	Color	Image	Width x Length	Thickness	Text
	Green		5W x 8.2L	2 Mil	
	Green		6W x 10L	2 Mil	
	Green		8W x 14L	2 Mil	

Medication Transfer Bags



- **FDA and USDA-compliant polyethylene:** Ensures safety and adherence to regulations.
- **Secure medication delivery:** Protects medication integrity with write-on labeling.
- **Tamper-evident closure:** Provides tamper-proof security and easy opening.

Ordering Information

Closure	Color	Image	Material	Width x Length	Thickness	Text
Tamper Evident Tape	Clear	White Block Imprint	Polyethylene	8W x 10L	2 Mil	Med. Trans. Bag

CUSTOM PACKAGING + SPECIALTY BAGS

Narcotic Return Bags



- **Safe and Secure:** Meets FDA and USDA standards for medical use. Tracks controlled substances throughout the dispensing and return process.
- **Efficient Documentation:** High-contrast design (black print on white) for easy record-keeping. Conveniently packaged in bulk (100 per case).

Ordering Information

Closure	Color	Image	Width x Length	Thickness	Text
Zip Seal Closure	Clear	Write-On Area	6.5W x 8L	2 Mil	

Patient Medication Bags



- **FDA and USDA-approved polyethylene:** Ensures patient safety and meets relevant regulatory standards.
- **Tamper evident bag:** Protects medication integrity with a clear indication of tampering.
- **Individually numbered:** Simplifies tracking and inventory management.

Ordering Information

Closure	Color	Image	Material	Width x Length	Thickness	Text
Tamper Evident Tape	Clear	White Block Imprint	Polyethylene	9W x 15L	3.25 Mil	Security Bag

CUSTOM PACKAGING + SPECIALTY BAGS

PRN Zip Bags




- **Certified for medical use:** Meets both FDA and USDA standards for safety and quality, ensuring your confidence in the product.
- **Clear medication identification:** Features orange printing to instantly alert staff of prescribed medication, promoting patient safety and reducing potential errors.
- **Bulk packaging for convenience:** Comes in cases of 100, offering efficient storage and minimizing the need for frequent restocking.

Ordering Information

Closure	Color	Image	Material	Width x Length	Thickness	Text
Zip Seal Closure	Clear		Polyethylene	3W x 5L	2 Mil	PRN
Zip Seal Closure	Clear		Polyethylene	4W x 6L	2 Mil	PRN
Zip Seal Closure	Clear		Polyethylene	6W x 9L	2 Mil	PRN

Refrigerate/Frozen Zip Bags



- **Medication temperature-sensitive:** Designed to clearly identify medications requiring refrigeration for safe storage and handling. Visual cues alert healthcare staff to ensure proper care.
- **Convenient bulk packaging:** Packed in cases of 100 for efficient inventory management and distribution.

Ordering Information

Closure	Color	Image	Material	Width x Length	Thickness	Text
Zip Seal Closure	Clear			4W x 6L	2 Mil	Refrigerate
Zip Seal Closure	Clear			6W x 9L	2 Mil	Refrigerate
Zip Seal Closure	Clear			9W x 12L	2 Mil	Refrigerate
Zip Seal Closure	Clear			12W x 15L	2 Mil	Refrigerate

CUSTOM PACKAGING + SPECIALTY BAGS

STAT Zip Bags



- **Instantly recognize critical moments:** A vibrant red imprint stands out, immediately alerting users to urgent situations.
- **Built for durability:** Double-walled construction with 2 Mil thickness ensures these supplies can withstand demanding environments.
- **Streamlined organization:** Inner packs of 100 keep your inventory tidy and accessible, ready for quick deployment.

Ordering Information

Closure	Color	Image	Material	Width x Length	Thickness	Text
Zip Seal Closure	Clear			4W x 6L	2 Mil	STAT
Zip Seal Closure	Clear			5W x 8L	2 Mil	STAT

Tamper Evident Bags



- **FDA and USDA-approved polyethylene:** Ensures patient safety and meets relevant regulatory standards.
- **Tamper evident bag:** Protects medication integrity with a clear indication of tampering.
- **Individually numbered:** Simplifies tracking and inventory management.

Ordering Information

Closure	Color	Image	Material	Width x Length	Thickness	Text
Tamper Evident Tape	Clear			6W x 6L	2 Mil	
Tamper Evident Tape	Clear			8W x 10L	2 Mil	
Tamper Evident Tape	Clear			10W x 14L	2 Mil	

CUSTOM PACKAGING + SPECIALTY BAGS

95kPa Specimen Transport Bags



- **Safeguard sensitive samples with the 95kPa PressurePro™ transport bag:** This leakproof specimen container boasts 95kPa pressure resistance, exceeding International Air Transport Association (IATA) standards for safe air travel of hazardous liquids. Its custom film construction and powerful adhesive closure create a continuous, impenetrable seal, ensuring critical liquid and substance containment throughout transport. Completely customizable.

Ordering Information

Closure	Color	Image	Material	Width x Length	Thickness	Text
Adhesive Seal Closure	Clear	Biohazard	Proprietary	8W x 10L	6 Mil	95kPa

CUSTOM PACKAGING + SPECIALTY BAGS

Custom Mailer Boxes

CUSTOM BOX SIZES

We produce your custom mailers in the exact size you need within the following ranges:

- Length: 3" - 25"
- Width: 2" - 25"
- Depth: 1" - 15"



The sizes we show correspond to the interior dimensions. Depending on your box's application, you may want to add some clearance to each side to make sure your products fit perfectly.

MATERIAL OPTIONS

To start customizing your printed corrugated box, choose from the following options:

Standard White Corrugated Cardboard

- Popular, economical option
- Made of sturdy, sustainable material
- Uncoated finish
- Improved print quality with HD Print

Premium White Corrugated Cardboard

- Smooth surface with a luxurious feel
- Clay-coated for a whiter brighter finish
- Best for luxury brands, gift boxes, and promotional kits
- Upgraded HD Print available exclusively with Satin Finish
- Select "Premium With Glossy Ink" for a high-gloss UV finish on printed areas

Kraft (Brown) Corrugated Cardboard

- Rustic brown kraft cardboard with a raw, natural look
- Ideal for brands promoting natural products & ingredients
- Greener single-pass HD Print process means less waste

HD Print

- Smaller ink dots bring out detail in your design
- Odorless, true water-based inks comply with even the most stringent global food safety regulations
- Aqueous coating gives the surface a consistent look and premium feel

MATERIAL THICKNESS

We automatically select a material thickness that works best for your box dimensions. If you require a specific thickness, please contact us.

- 1/16" E-Flute Corrugated Cardboard.
- 1/8" B-Flute Corrugated Cardboard.


CUSTOM PACKAGING + SPECIALTY BAGS

Custom Product Boxes

CUSTOM BOX SIZES

We produce your custom mailers in the exact size you need within the following ranges:

- Length: 0.75" – 30"
- Width: 0.75" – 20"
- Depth: 1.75" – 20"



Note that these measurements refer to the interior or inside area of the box. You may want to add some clearance to each side of your product dimensions, depending on how you want it to fit inside the box.

MATERIAL OPTIONS

Cardstock

This durable material is thicker than paper but thinner than cardboard. The larger the stock number, the thicker the stock.

- Offers the highest print quality for in-store displays
- Your choice between 3 different coatings for the feel and look you want
- 18 pt. (thickest cardstock) is recommended for larger boxes

Corrugated Cardboard (1/6" E-Flute)

This is a thicker material that creates a more durable product box.

- Stands up to rough handling
- Can go directly in the mail if needed
- White stocks have a soft, smooth surface with a luxurious feel

PREMIUM FEATURES

Add a special touch to make your packaging extraordinary.

- Embossing
- Spot UV
- Foil stamping
- Soft touch

To add a special feature to your boxes, contact us for a quote.


CUSTOM PACKAGING + SPECIALTY BAGS

Custom Shipping Boxes

CUSTOM BOX SIZES

Your order will be produced to the exact size you need within the following ranges:

- Length: 2" – 30"
- Width: 2" – 29"
- Depth: 2.5" – 32"



Note that the customizable measurements refer to the interior dimensions of the box. You may want to add space on either side of your custom retail packaging when determining the size of your box to leave room for packaging or inserts.

MATERIAL OPTIONS

Standard White Corrugated Cardboard

- Popular, economical option
- Made of sturdy, sustainable material
- Uncoated finish
- Improved print quality with HD Print

Premium White Corrugated Cardboard

- Smooth surface with a luxurious feel
- Clay-coated for a whiter brighter finish
- Best for luxury brands, gift boxes, and promotional kits
- Upgraded HD Print available exclusively with Satin Finish
- Select "Premium With Glossy Ink" for a high-gloss UV finish on printed areas

Kraft (Brown) Corrugated Cardboard

- Rustic brown kraft cardboard with a raw, natural look
- Ideal for brands promoting natural products & ingredients
- Greener single-pass HD Print process means less waste

HD Print

- Smaller ink dots bring out detail in your design
- Odorless, true water-based inks comply with even the most stringent global food safety regulations
- Aqueous coating gives the surface a consistent look and premium feel

MATERIAL THICKNESS

We automatically select a material thickness that works best for your box dimensions. If you require a specific thickness, please contact us.

- 1/16" E-Flute Corrugated Cardboard.
- 1/8" B-Flute Corrugated Cardboard.



Cenmed Enterprises, Inc.

121 Jersey Avenue | New Brunswick, NJ 08901 | USA

Customer service, sales and technical support:

Phone | (732) 447-1100

Email | info@cenmed.com

Website | www.cenmed.com